





Flexible and Wearable sensors for health monitoring applications

Reva Teotia (2019A3PS0268P)

Introduction


Wearable devices for continuous and remote health monitoring.

Requirements

- Light weight
- Flexible
- Biocompatible

Applications

- Pulse rate measurement
- Motion detection
- Glucose, pH and other biological indicators detection
- Physiological signal like EEG, EMG, ECG monitoring




Temperature Sensors

Temperature sensors are vital for monitoring the body temperature. Temperature outside the normal range can be critical for the health.

Transduction mechanism

- Flexible resistance temperature detector (FRTD)
- Thermistors
- Thermocouple



innovate

Temperature Sensors

Flexible resistive temperature detector (FRTD)

Change in temperature cause a change in electrical resistance.

- Temperature coefficient of resistance (TCR) = $(\Delta R/R_0)/\Delta T$
- Sensitivity $S = \Delta R/R_0$

Active material - Metals like Ni, Cu, Pt

Carbon derivatives - CNT,


graphene, quantum dots

Conducting polymers -

PEDOT:PSS

Fabrication methods -

Laser Digital Printing (LDP) Inkjet printing Sputtering Spin coating



lead

Thermistors


The resistance changes occurs due to change in temperature.

Active material - Metals like Pt, Cu, Ag, Au, AgNW Carbon derivatives like rGO, CNT, graphene

Substrate - PI, PDMS, PET, PEN

Fabrication methods -

Inkjet printing Sputtering Spin coating



Thermocouples


Seebeck effect - production of emf in two dissimilar conductors when the junctions are maintained at different temperatures.

Active material - alloy films

- nickel-aluminum-silicon-manganese alloy film,
- nickel-aluminum alloy film,
- p-Sb2Te3 film, n-Bi3Te3 film, Bi-Te film and Sb-Te film

Fabrication methods -

Sputtering Spin coating



Heat Sink

K-type Thermocouple

Heater K-type

Thermocouple

+Heater

Thermal pad

+Heat Sink


innovate

Pressure Sensors

Pressure sensors can be used for determining pulse, blood pressure or external pressure applied on e-skin or prosthetics.

Transduction mechanism

- Capacitive
- Piezoelectric
- Piezoresistive




Capacitive Pressure sensor

$$C = \varepsilon_r \varepsilon_0 A/d$$

Where ε_r is the relative permittivity, A is the plate area and d is the gap between plates.

Active material - dielectric material - PDMS, PU, organosilicon elastomer, polystyrene, PET, woven-spacer structures


Microstructures are created to enhance sensitivity.




Pressure Sensors

Piezoelectric Pressure sensor

Piezoelectric effect - The application of external pressure that cause deformation generates current and vice versa.






Active material - P(VDF-TrFE)), barium titanate (BaTiO3), lead zirconate-titanate (PZT), and zinc oxide (ZnO)


Pressure Sensors

Piezoresistive Pressure sensor

Piezoresistance - Deformation in the structure cause change in electrical resistance.

Active material - metallic particles, metal NW carbon-based materials including carbon black, CNTs, and reduced graphene oxide (rGO)





Strain Sensors

Strain is the relative change in the length of a structure under stress.

Strain sensors convert mechanical stimuli into optical or electrical signals.

Application - Recognize the posture or gesture

Continuous blood pressure monitoring

Artificial tactile sensation


Vocal-cord vibrations


Pulse


Respiration


Transduction mechanism

- Piezoresistive
- Capacitive
- Piezoelectric










Piezoresistive Strain sensor

Tensile strain Gauge factor (GF) = $(\Delta R / R0) / \epsilon$, that is, the fractional change in the resistance, with ϵ being the strain



Piezoresistive strain sensor

Compressive strain Sensitivity = $(\Delta E / E0) / \Delta P$. where ΔE represents the change in the electrical signals

Active material - CNT, graphene, carbon black, polymers, metal nanowires

Advantages - simple structure, low cost, high linearity and excellent sensing performance




Strain Sensors

Capacitive Strain sensor

Gauge factor, $GF = (\Delta C / C0) / \epsilon$

Active material - dielectric material - elastomers, ionic liquid/gels, thinned inorganic material

Advantages - Low input energy, High dynamic response and Low noise impact



innovate

Piezoelectric Strain sensor

These can be used as a generator as well as the biosensor.

Evaluated in either of the two working modes - open circuit voltage (V_{oc}) Short circuit current (I_{sc})

Active material - ZnO, PZT nanoribbons, BaTiO3, and organic polymer like PVDF

Thank You